



FAMO100 Moisture Analyzer Touch Screen

External Calibration



- Capacity: 100g
- Readability - 1mg
- Repeatability - ± 0.2mg
- Linearity - ± 0.2mg
- Pan Size - Φ100mm
- Output Interface - RS232C
- Dimensions (L * W * H) - 340mm × 215mm × 350mm
- Net Weight - 6000g
- Power Supply - Universal Adaptor AC 110-240V.
- Display: LCD with Backlight, Touch Screen
- Unit: Moisture, Dry,

FAMO Fast and accuracy moisture determination. Dry weight range determination from 100% to 0



External Calibration

External Calibration Calibration Weight Included



Touch Screen Display

Capacitive Touch Screen, High Resolution,



Software for data collection

software for communication with the scale provides the ability to collect data, control the measurement process (optional)

TYPE	FAMO100
PRICE	
Calibration	External
Measuring range g	100
Temperature Range	60°C-180°C
Resolution	1mg
Pan size	100mm
Display	High Resolution true color 5,7 inch LCD
System	Moisture Analyzer 100,00%-0.00%