



TECB Industrial Precision Balance with Remote Display

External Calibration

Download
Catalogue



TECB Series High Loading Balances have excellent weighing performance and perform routine and research weighing. With all range custom calibration method, weighing in all range is sharp and stable. Backlight increases visibility in all conditions. Optional battery power provides the opportunity to work with scales independent of the power supply.

- Capacity: 6/10/15/20kg/30kg
- Resolution: 0.1g
- Pan Size: 310 x 220mm
- Display: LCD with Backlight

TECB High Capacity Precision Balance with Detachable Display



External Calibration

External Linear Calibration



Weighing - Piece Counting

Different type of weighing mode,



Detachable Display

Easy manipulation, Beck light display ...

30kg/1g

**External
Calibration**

Big Pan Size

310*220mm

**Detachable
Display**

Piece Counting

Type	TECB6	TECB10	TECB15	TECB20	TECB30
Maximum Capacity (kg)	6	10	15	20	30
Accuracy Class	III				
Division (g)	1g				
Stabilization time	≤2sec				
Operating temperature	15°C~30°C, fluctuation range <5°C/h				
Operating Relative humidity (RH)	40%~80%				
Power supply	AC110-240V 50/60Hz DC5V-2A (adaptor)				
Pan Size (mm)	255*190mm				
Calibration	External				
Display	LCD with backlight				
Interface	RS232C (Optional)				
System	Check Weighing, Parts Counting, Percentage Weighing				
Optional Accessories	Rechargeable battery, Under Hook, Double Display, Printer, Instruction Manual,				
Dimension(L*W*H)	32×26×11 (L*W*H)				
Gross weight kg	4				

Please note

If this Balance does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts

Laboratory Equipment production and Distribution

www.colo.si

Presernov Trg 9

8000 Novo Mesto

Slovenia

tel. +386 590 48 880 E-mail: office@colo.si

www.colo.si

© 2021 COLO Distributors Price List V1