



LACE12RM Refrigerated High Speed High Capacity Centrifuge



LED Display



Refrigerated
-20°C to 40°C

Download
Catalogue



Angle Rotor

- 30122 Angle Rotor 6x300ml 12000rpm/22650xg
- 30123 Angle Rotor 6x500ml 8000rpm/11740xg
- 33097 Angle Rotor 4x1000ml 4000rpm/4060xg



Max. Capacity 4x1000ml

Angle Rotor



12000rpm/19040xg

Maximum speed and rcf meet all demand for scientific, clinical PCR centrifuge,



Angle rotors

Adapters available for this rotor



LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

Item Name	LACE-12RM
Max. Rotor capacity	4x1000ml
Max. Speed rpm	12,000rpm/22,650xg
Noise	<58dB (A)
Timer	99H, 59min
Rotor type	Angle rotor
Power supply	AC 220V/50Hz
Product dimension	810 x 650 x 870 mm
Weight kg NV/GW	235kg/275kg