



# LACE16RH Refrigerated High Speed Centrifuge

## Angle Rotor



LED Display



Refrigerated  
-20°C to 40°C

Download  
Catalogue



### Max. Capacity 6x100ml

Angle Rotor



### 16000rpm/20600xg

Maximum speed and rcf meet all demand for scientific, clinical PCR centrifuge,



### Angle rotors

Adapters available for this rotor



### LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



### Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



### Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

- 30402 Angle Rotor 40x0.5ml 14000rpm/19970xg
- 30401 Angle Rotor 12x1.5/2ml 16000rpm/16260xg
- 30403 Angle Rotor 24x1.5/2ml 15000rpm/20600xg
- 30404 Angle Rotor 30x1.5/2ml 13500rpm/19340xg
- 30405 Angle Rotor 16x5ml 15000rpm/19350xg
- 30406 Angle Rotor 12x7ml 14000rpm/16370xg
- 30407 Angle Rotor 12x15ml 10000rpm/11840xg
- 30408 Angle Rotor 12x10ml 13000rpm/14510xg
- 30409 Angle Rotor 8x20ml 12000rpm/14510xg
- 30410 Angle Rotor 6x30ml 12000rpm/14000xg
- 30411 Angle Rotor 6x50ml 11000rpm/13480xg
- 30412 Angle Rotor 6x70ml 10000rpm/10810xg
- 30413 Angle Rotor 4x100ml 12000rpm/14850xg
- 30414 Angle Rotor 6x100ml 8000rpm/7280xg
- 30416 Angle Rotor 30x0.5ml 15000rpm/18510xg
- 30444 Angle Rotor 48x1.5/2ml 11000rpm/12840xg