





# LACE24M Universal Laboratory Centrifuge with angle rotor 24x10ml

-  LED Display
-  Neutral Temperature

[Download Catalogue](#) 

## Angle Rotor

Laboratory centrifuge LACE24 with the 24x10ml angle rotor, is an economical solution for smaller laboratory work. Small size and weight allow for easy transportation from one place to another. LACE24 Centrifuge is ideal as an additional instrument in a laboratory, serving either as a reserve or for sample preparation.

## LACE24M



**4000rpm/2250xg**  
Maximum speed and rcf meet all demand for clinical centrifuge



**Angle Rotor 24x10ml**  
Combined with adapter this type of rotor can work with all type of centrifuges from 2ml to 15 ml,



**Digital Display**  
Speed/RCF, Time, Parameters can be set during the centrifuge work..



**Brushless Motor**  
COLO brushless motor is maintenance free, no carbon dust pollution



**Safety Function**  
Electric Lead Lock with safety function, disbalance detection, over heating of centrifuge protection

- 31602 - Angle, 12x15ml - 4000rpm/2150xg
- 31604 - Angle, 18x10ml - 4000rpm/2250xg
- 31607 - Angle, 24x10ml - 4000rpm/2220xg
- 31601 - Angle, 30x5ml - 4000rpm/2010xg
- 32665 - Swing out, 4x10ml - 4000rpm/2610xg
- 31606 Angle Rotor 6X50ml 4000rpm 2100xg

330€

330€