





LACE32 Low Speed Centrifuge with 32x10ml Swing out Rotor

-  LED Display
-  Neutral Temperature

Download Catalogue 

Angle/Swing out Rotor

- 30569 - Swing out, 4x250ml - 4000rpm/2880xg
- 31376 - Swing out, 24x15ml - 4000rpm/2950xg
- 31377 - Swing out, 32x10ml - 4000rpm/2880xg
- 31378 - Swing out, 48x5 ml - 4000rpm/2480xg
- 31379 - Swing out, 48x7ml - 4000rpm/2760xg
- 30614 - Angle, 6x100mk - 5000rpm/3130xg
- 31380 - Swing Micro plate, 2x2x96well -
- 30641 - Angle, 12x50 - 5000rpm/3860xg
- 30639 - Angle, 24x15ml - 5000rpm/3500xg



Max. Capacity 4x250ml

Swing out Rotor



5000rpm/4807xg

Maximum speed and rcf meet all demand for clinical centrifuge



Swing our Rotor 32x10ml/15ml

Is included as a standard configuration,



TOUCH Screen Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust polution



Safety Function

Electric Lead Lock with safety function, disbalance detection, over heating of centrifuge protection

Item Name	LACE-32
Max. Rotor capacity	4X250ml
Max. Speed rpm/xg	5000rpm/4807xg
Power	450W
Ambient temperature	0-40°C
Noise	≤55dB (A)
Timer	0-99h59nin / continuous / short spin button (option)
Electric safety lid lock	Mechatronic motor door lock
RPM/RCF conversion	Yes
Acceleration/Brake	10 programming step of acceleration and deceleration
Power supply	AC 220V/50Hz
Weight kg NV/GW	37/41