



LACE-6AH Low Speed Universal Centrifuge with rotor 6x50ml



LED Display



Neutral
Temperature

Download
Catalogue



Angle/Swing out Rotor

LACE-6AH is High Capacity, economical universal centrifuge suitable for large clinical laboratories.

Included



Max. Capacity 4x180ml

Swing out Rotor



6000rpm/4430xg

Maximum speed and rcf meet all demand for clinical centrifuge



Angle rotor 6x50ml

Adapters available for this rotor



LED Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



Safety Function

Electric Lead Lock with alarm with safety function, disbalance detection, over heating of centrifuge protection

Item Name	LACE-6AH
Max. Speed	6000rpm/4110xg
Capacity (Angle rotor)	6 x 50ml 12x15ml Adapters conical tube 17,5x125 Vacuum tube 10ml 17x107mm
Capacity (Swing out rotor)	4 x 180ml
Noise	<45dBA
Timer	0-99 min
Rotor type	Swing out and Angle Rotor
Power supply	AC 220V
Product dimension	540 x 370 x 280mm
Weight kg NV/GW	33kg/37kg