



MICE8 Small Laboratory Centrifuge 4000rpm/2800xg

Little Giant



LED Display



Neutral Temperature

Download Catalogue



Laboratory centrifuge MICE8 with the eight place fixed rotor, is an economical solution for smaller laboratory work. Small size and weight allow for easy transportation from one place to another. MICE8 Centrifuge is ideal as an additional instrument in a laboratory, serving either as a reserve or for sample preparation.

MICE8 Small Laboratory Centrifuge light weight, small size, proper speed for blood , urine PRP separation...



4000rpm/2800xg

Maximum speed and rcf meet all demand for clinical centrifuge



Angle Rotor 8x15ml

Combined with adapter this type of rotor can work with all type of centrifuges from 2ml to 15 ml,



Digital Display

Speed/RCF, Time, Parameters can be set during the centrifuge work..



Brushless Motor

COLO brushless motor is maintenance free, no carbon dust pollution



Lead Lock

Simple Lead Lock with one button

PRODUCT CODE	MICE8
PRICE	
Speed	4000rpm/1800xg
Capacity	8x15ml
Timer	0-99min
Noise	<45dBA
Display	LED - Display: Speed, RCF, Time
Power supply	220V (±10%) 50-60Hz
Safety	Mechanical lid lock
Weight NW/GW kg	7/8,5kg
Product Dimension	260×265×310mm