



LACETRAP Cold Trap for Vacuum Concentrator

LACETRAP

Download
Catalogue



Low Noise



Air colled



Tank capacity 7 L



Fast cooling



Aggressive resistant

Suitable for highly aggressive or corrosive gases and vapor (series COT and CDT)

Function	Vacuum
Tank capacity	7L
Min. Temperature range	≤35°C
Temperature sensor	PT100 class A
Cooling system	CFC compressor cooled by air
System control	Self-setting function, Microprocessor digital control on/off START/STOP Timer
Timer range	Continues, 99min/59 sec, 99h/59min, 99 days 23h
Chamber	seamless stainless steel 304
Accuracy	±2°C
Temp resolution	1°C
Relative humidity	(60-80%)Ventilation
Cooling Coil	Copper(Nickle Plating)
Power supply	220V±10% 50HZ
Dimension mm	370x460x640
GW	53kg

Please note

If this Vacuum Pump does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts



Laboratory Equipment production and Distribution

www.colo.si

Presernov Trg 9, 8000 Novo Mesto, Slovenia

tel. +386 590 48 880 E-mail: office@colo.si

www.colo.si

© 2021 COLO Distributors Price List V1