



ULT950 -40 to -86C ULT Freezer 500L

Double layer door structure

LCD



Smart air
circulation

Download
Catalogue



Interior made from ASI304 stainless steel
Noise level less then 52dBa

-4°C to -86°C



LCD Touch Screen

Microcomputer control temperature and oscillation frequency with timing functions, Built-in power-off protection can automatically resume operation after power on.



USB port to upload and download data



Password protected data



Optional cloud temperature control



Optional electronic lock data card

- Ultra Low Temperature Freezer
- Features
- Temperature range: -40~-86°C
- 500 L Storage capacity
- HC refrigerant, environmental friendly
- Microprocessor controller
- LCD touch screen for intuitive operation
- SUS/AISI304 interior working chamber, epoxy coated exterior, PURF insulation
- Test hole QTY/diameter 1/25MM
- Shelf/drawer QTY : Stainless steel shelf/3
- Full loading: racks quantity 16(5X5), Boxes (2 Inches) 400pcs, Samples(2ml) 40000
- Password protect the freezer from unauthorized access
- 72 hours power failure backup system (alarm)
- 891×1155×1994.5 Exterior size(W*D*H)(mm)
- 307/368 NT./GT. (kg)

Remote alarm device and system for ULT freezer

Feature

Temperature range: -200~150°C

Temperature accuracy: ±0.5°C

Data recording interval: 1min-24hours

Data recording rapacity: 100000

Data transmitting way: 4G/WIFI/local server

Battery: rechargeable

Optional accessories:

Aluminum racks, box, electornic loc cards,

Please note

If this Ifreeezer does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

www.colo.si

© 2023 COLO Distributors Price List V3

COLO LabExperts

Laboratory Equipment production and Distribution

www.colo.si

Polje ob Sotli 4, 3250 Buce

Slovenia

tel. +386 64 222 724 E-mail: office@colo.si

