



GOBI – Laboratory Hotplate Aluminium Surface

 Max. temp.
350°C

Download
Catalogue



t° LIMITER

- COLO Heating plate is ideal for precise temperature control when heating media. The three base models provide all the capabilities your lab needs for heating. Short heating times, precise temperature control, chemical resistance are just some of the benefits of COLO Hot Plates.
- Chemical resistant coating
- Digital model equipped with PID control to prevent overheating of medium
- Auto tuning function for maximum performance of PID control and heating time
- All models are equipped with heat indicator
- A separate On/Off button for the heating function prevents unintentional heating



Temperature range RT+5°C to 350°C

High temperature uniformity,



Intelligent PID Temperature Control

Digital microprocessor PID controller with integrated auto diagnostic system and error indicator



Special coating aluminium surface

Aluminium surface with perfect temperature uniformity,



LED Display to read and setting

Simple and easy to navigate, read and import parameters



Timing Function

Timer with time delay function

Hotplate type	GOBI5A	GOBI8A	GOBI9A	GOBI5D	GOBI8D	GOBI9D	THAR5D	THAR8D	THAR9D	
Price										
Controller Type	Analogue			Digital microprocessor PID control, Correction of temperature probe deviation, PID temperature control, high temperature limit, timer,						
Heating surface	Aluminium						Graphite			
Heating surface dimension mm	300*300	300*450	400*600	300*300	300*450	400*600	300*300	300*450	400*600	
Maximum temperature	320°C	320°C	320°C	350°C	350°C	350°C	450°C	450°C	450°C	
Heating power	1200W	1800W	3000	1200W	1800W	3000	1200W	1800W	3000	
Temperature sensor	Thermostat			K type A class						
Temperature accuracy	±1%							±10%		
Plate temperature uniformity	±10%							±15%		
External size	34*37*22	42*50*22	47*65*22	34*37*22	42*50*22	47*65*22	34*37*22	42*50*22	47*65*22	
Plate load	20kg	30kg	40kg	20kg	30kg	40kg	20kg	30kg	40kg	
Timer	Digital							X		
Voltage	220/230V 50-60Hz									
NW/GW kg	6,2/7	8,3/9,6	12/14,3	6,2/7	8,3/9,6	12/14,3	6,2/7	8,3/9,6	12/14,3	

Please note
 If this Hotplate does not match the specific needs of your application, or some options are not listed for sale, please feel free to contact us. Our manufacturing engineers will come up with technical solutions to meet your needs. We reserve the right to change technical specifications at any time.

COLO LabExperts 
 Laboratory Equipment production and Distribution
 www.colo.si
 Presernov Trg 9, 8000 Novo Mesto, Slovenia
 tel. +386 590 48 880 E-mail: office@colo.si