

ORSHA410 Rocker Shaker 10-70rpm

Product Details


MANUFACTURER:

COLO CODE: PRO5323

PRODUCT CATEGORY: Shakers

PRODUCT SUBCATEGORY:

TYPE: Rocker



Technical Specification

Features:

- Double LCDs for independent display of time and speed.
- Different platforms and fixtures for selection.
- 3D shaker: gentle 3D motion with fixed angle 9°; Rocking shaker: smooth see-saw rocking motion with fixed angle 9°.
- Wide range of platforms for use with a variety of vessels.

Technical Parameters:

Model	PRO5323	PRO4883
Item No.	8032422100	8032322100
Rocking Motion	3D orbital	See-saw rocking
Angle	0-15°	
Max. Load Capacity	5kg	10kg
Motor Type	DC motor	
Speed Range	10~70rpm	
Speed And Timer Display	LCD	
Time Range	1~1199min	
Operation Mode	Continuous/ timed operation	
Working Environment	5-40°C, 80%RH	
Motor Input	IP21	
Motor Output	40W	
Protection Type	24W	
Consumption	50W	
Power Supply	100~240V,50/60Hz	
Standard Accessory	1 pc SK330.5 petri dish platform	
External Size(W*D*H)	360*430*106mm	360*410*200mm
Package Size(W*D*H)	540*480*370mm	540*480*340mm
Gross Weight	16kg	15kg

Optional Accessories