

ORSHA510 Rocker Shaker 1-80rpm

Product Details

Product Number

PRO5821



Technical Specification:

Technical Parameters :

Model	ORDHA530
SKU	PRO5821
Max load	5KG
Motor type	DC brushless motor
Motor input power	35W
Motor Output power	25W
Speed range	1~80rpm
Slope	10°
Speed display	LCD
Timer display	LCD
Time set range	99h59min
Power supply	100~240V,50/60HZ
NW	8.5KG
Allowable ambient temperature/humidity	5~40°C / 80%
Dimension L*W*H mm	330*370*140
Shell Protection Class	IP21
Data	RS232 optional

Accessories:

Universal Tray Size:210*282mm

Adjustable tray Size:210*282mm

Spring tray Size:210*282mm

Fixed tray Size:210*282mm

Flask tray 6*250ml

Flask tray 8*100

Flask tray 8*50

Flask tray 8*25

Fixed tray (Double layer) Each layer size:210*282mm

Universal tray (Double layer) Each layer size:210*282mm

Gil pad Size: 205*275mm, used with 210*282mm universal tray

Custom tray