



VOMI25D High Speed Laboratory Vortex Mixer Digital Control

Price - €

Download
Catalogue



- The High Speed laboratory Vortex Mixer provides smooth instant vortex of even the largest sample sizes.
- Available for different turn plate and laboratory containers.
- Functions of vortex mixing - orbital mixing
- Orbital diameter 4 mm.
- Digital timer 1-999min
- Choice of continuous rotating and touch rotating mode
- Triple mode switcher – on /off/touch
- Multi function Vortex Mixer
- Variable speed Set, 0 - 4200rpm,
- Pop-off Cup Head for test tubes up to 50ml (30mm)
- Shaking Motion: Orbital
- Control: Digital knob adjustable
- Dimensions (mm): 138×136×205
- Weight: 2,8 kg
- Power: AC 100-240 V, 50 / 60 Hz - 15W . (±10%)



Digital Speed Regulation

0-4200 rpm



Continuous /Touch rotating operation



DC brushless motor

Maintenance free



Orbital shaking motion 4mm



Digital Display