



VOMI25S High Speed Laboratory Vortex Mixer

Price - €

Download
Catalogue



- The High Speed laboratory Vortex Mixer provides smooth instant vortex of even the largest sample sizes.
- Available for different turn plate and laboratory containers.
- Functions of vortex mixing - orbital mixing
- Orbital diameter 4,5 mm. $\pm 10\%$
- Choice of continuous rotating and touch rotating mode
- Rotate speed by electronic adjustable
- Multi function Vortex Mixer
- Variable speed Set, 0 - 4500rpm,
- Pop-off Cup Head for test tubes up to 50ml (30mm)
- Shaking Motion: Orbital
- Control: Analogue knob adjustable
- Dimensions (WxDxH) (mm): 127x130x160
- Weight: 0,8 kg
- Power: AC 100-220 V, 50 / 60 Hz - 15W . ($\pm 10\%$)



Analogue Speed Regulation

0-4500 rpm



Continuous /Touch rotating operation



DC brushless motor

Maintenance free



Orbital shaking motion 4,5mm $\pm 10\%$



Compact design