



VOMI55 Orbital Micro-Plate Mixer

Price - €

Download
Catalogue



Analogue Speed Regulation

0-1500 rpm



Continuous operation



DC brushless motor

Maintenance free



Wide range of accessories



Compact design

Item Name	VOMI55
Shaking Mode	Orbital
Max. Speed rpm	0-1500
Orbital Diameter	4,5mm
Max shaking weight	0,5kg
Motor type	DC brushless motor
Motor input	18W
Motor Output	10W
Ambient working temp.	5-40°C
Relative Humidity	80%
Protection Type	IP21
Power consumption	20W
Power supply	C110/220V±10%, 50/60Hz
External size W*D*H mm	150*260*80
Gross weight kg	3