



UCB7 Ultrasonic Cleaning Bath with Heating

Price - €

Download
Catalogue



Digital Display

...



Timer

Timing range 1-99 min adjustable



Stainless Steel Basket

Designed especially for ultrasonic cleaner to reduce ultrasound cover



Stainless Steel Structure

The shell and cover, made from stainless steel isolated to decrease ultrasound noises



Heating

Chamber volume	6.5L
Number of transducer	3
Timer	1-99min
Heating power	300W
Capacitor (pF)	3500-3900 +/- 10%
Chamber thickness	0.8mm
Temperature range	Up to 80°C
Ultrasound frequency	40kHz
Transducer power	60W
Ultrasound power	180W
Basket dimension W*D*H mm	270*145*110
Water tank dimension W*D*H mm	300*150*150
Cleaner dimension W*D*H mm	330*175*320
Discharge valve	Yes
PID temperature regulation	Yes, up to 80°C
Auto stop	Yes
Power supply (euro plug)	220V/50Hz